

# HOW TO SELECT A POWER METER

By Burt Mooney  
Ophir-Spiricon



# How to select a power meter

- The first thing one generally has to do is look in depth at what type of laser has, and its parameters that one needs to measure.
- From that information you can determine what type of sensor (detector) one needs and then look at the instrument one will need to display the results.

# Determining Sensor Type

- Is laser Continuous Wave (CW)?
- Is laser Pulsed?
- Wavelength?
- Min and Max Energy or Power?
- Beam size

# Determining Sensor Type

- If the laser is continuous you will either need a thermopile or photodiode sensor.
- Knowing the min and max power and wavelength range will determine which of these to choose.



# Determining Sensor Type

- Thermopiles
  - Spectrally broadband from UV to FAR IR
  - 50 microW to 10 kiloW
- Photodiodes
  - Not spectrally flat, limited spectral range
  - Silicon 190 -1100 nm
  - Germanium 700 – 1800 nm
  - InGaAs 800 – 1700
  - picoW – 1 W

# Determining Sensor Type if Thermopile

## ➤ **Two Types of Absorbers**

### ➤ **Surface Absorbers**

- The laser light is absorbed in the front surface of the sensor.
- Optical black paint: 500 W/cm<sup>2</sup>, 50 mJ/cm<sup>2</sup> @ 10 nsec
- High Temp Ceramic: 26 kW/cm<sup>2</sup>, 300 mJ/cm<sup>2</sup> @ 10 nsec

### ➤ **Volume Absorbers**

- The laser light is absorbed in a bulk material and then conducted to the metal disc.
- Typical lasers include Q-switched relative high energy Nd:YAG, Ruby, Alexandrite, and

# If Your Laser is Pulsed

- You may need a pyroelectric sensor if:
  - Rep rate is single shot to 25 kHz
  - Spectrally broadband from UV to FAR IR

# Watt's a Joule?

- Conversely, if one knows the Average Power and repetition frequency, one can get the Average Energy per pulse.
- Energy = Average Power divided by pulse repetition frequency, i.e.  $1 \text{ mJ} = 1 \text{ W} / 1 \text{ kHz}$ .





# Watt's a Joule?

- **Peak Power**
  - Watts of Peak Power = Energy in Joules divided by pulsewidth (FWHM), For example:
  - 1 Joule @ 10 nsec = 100 MegaWatts
- **Fluence**, or energy density per unit area, is just Energy divided by area being irradiated.  
If energy is known just divide by area.  
For example: 1 mJ pulse in a 1 cm<sup>2</sup> area = 1 mJ/cm<sup>2</sup>
- **Irradiance**, or power density per unit area, is just Power divided by area

# Selecting an Instrument Readout

- What features do you need?
  - Statistics on board?
  - Portability?
  - Battery Operation?
  - Computer interface?
  - User interface?
  - Multi-channel?
  - Wireless?

That's It